

DNAJC6/AUXILIN antibody



Catalog Number: 109951

Product name

DNAJC6/AUXILIN antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:20-1:200

Antibody description

DNAJC6/AUXILIN Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human brain tissue, mouse brain tissue, rat brain tissue. Positive IHC detected in human brain tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 100-120kDa, 86kDa

Validations




Preparation

This antibody was obtained by immunization of DNAJC6/AUXILIN recombinant protein (Accession Number: NM_001256864). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109951(DNAJC6 antibody) at dilution of 1:800

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:109951(DNAJC6 Antibody) at dilution of 1:50 (under 10x lens)

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:




DNAJC6/AUXILIN antibody



Catalog Number: 109951

Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:109951(DNAJC6 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:109951(DNAJC6 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)