

DNAJC9 antibody



Catalog Number: 109953

Product name

DNAJC9 antibody

Specificity

Human; other species not tested.

Antibody description

DNAJC9 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, Raji cells, SH-SY5Y cells. Observed molecular weight by Western-blot: 30-35 kDa

Preparation

This antibody was obtained by immunization of DNAJC9 recombinant protein (Accession Number: NM_015190). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109953(DNAJC9 Antibody) at dilution of 1:1000