

DNASE2 antibody

Catalog Number: 109960

Product name

DNASE2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DNASE2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of DNASE2 recombinant protein (Accession Number: NM_001375). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109960(DNASE2 antibody) at dilution of 1:300



IP Result of anti-DNASE2 (IP:Catalog No:109960, 4ug; Detection:Catalog No:109960 1:300) with HeLa cells lysate 3200ug.