

DNM2 antibody



Catalog Number: 109963

Product name

DNM2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

DNM2 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, HeLa cells, human brain tissue, Jurkat cells, mouse brain tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 100kd

Preparation

This antibody was obtained by immunization of DNM2 recombinant protein (Accession Number: NM_004945). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:109963(DNM2 antibody) at dilution of 1:500