

DNMT2 antibody



Catalog Number: 109967

Product name

DNMT2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DNMT2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, MCF7 cells, mouse liver tissue. Observed molecular weight by Western-blot: 43-45 kDa

Preparation

This antibody was obtained by immunization of DNMT2 recombinant protein (Accession Number: NM_004412). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109967(DNMT2 antibody) at dilution of 1:600