

# DNMT3A antibody



Catalog Number: 109968

## Product name

ELISA, IHC

DNMT3A antibody

## Dilutions

## Specificity

Recommended Dilution:

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200


## Antibody description

DNMT3A Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue, human kidney tissue.

## Validations

## Preparation

This antibody was obtained by immunization of DNMT3A recombinant protein (Accession Number: XM\_017003526). Purification method: Antigen affinity purified.




## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Immunohistochemical of paraffin-embedded human heart using Catalog No:109968(DNMT3A antibody) at dilution of 1:50 (under 10x lens)

## Storage

Store at -20°C. DO NOT ALIQUOT



## Clonality

Polyclonal

## Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human heart using Catalog No:109968(DNMT3A antibody) at dilution of 1:50 (under 40x lens)

## Applications