

DNMT3A antibody



Catalog Number: 109968

Product name

DNMT3A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DNMT3A Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue, human kidney tissue.

Preparation

This antibody was obtained by immunization of DNMT3A recombinant protein (Accession Number: XM_017003526). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human heart using Catalog No:109968(DNMT3A antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:109968(DNMT3A antibody) at dilution of 1:50 (under 40x lens)