

DNMT3L antibody



Catalog Number: 109970

Product name

DNMT3L antibody

Specificity

human; other species not tested.

Antibody description

DNMT3L Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, human testis tissue. Observed molecular weight by Western-blot: 44kd and 70kd

Preparation

This antibody was obtained by immunization of DNMT3L recombinant protein (Accession Number: NM_175867). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109970(DNMT3L antibody) at dilution of 1:300