

DNTT,TDT antibody

Catalog Number: 109972

Product name

DNTT,TDT antibody

IP: 1:200-1:1000

IHC: 1:20-1:200


Immunogen

[Human DNTT,TDT Recombinant protein \(His tag\)](#)

Validations

Specificity

Human; other species not tested.




Antibody description

DNTT,TDT Rabbit Polyclonal antibody. Positive IP detected in Raji cells. Positive IHC detected in human lymphoma tissue.

Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:109972(DNTT,TDT antibody) at dilution of 1:50 (under 10x lens)

Preparation

This antibody was obtained by immunization of DNTT,TDT recombinant protein (Accession Number: NM_004088). Purification method: Antigen affinity purified.



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:109972(DNTT,TDT antibody) at dilution of 1:50 (under 40x lens)

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, IHC, IP

Dilutions

Recommended Dilution:



DNTT,TDT antibody



Catalog Number: 109972

IP Result of anti-DNTT,TDT (IP:Catalog

No:109972, 4ug; Detection:Catalog No:109972
1:300) with Raji cells lysate 2000ug.