

DISP1 antibody



Catalog Number: 109979

Product name

DISP1 antibody

Specificity

Human; other species not tested.

Antibody description

DISP1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 150-171 kDa

Preparation

This antibody was obtained by immunization of DISP1 recombinant protein (Accession Number: BC011542). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109979(DISP1 antibody) at dilution of 1:500