

DKC1 antibody



Catalog Number: 109981

Product name

DKC1 antibody

Specificity

Human; other species not tested.

Antibody description

DKC1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, COLO 320 cells, HeLa cells, Jurkat cells, MCF-7 cells. Observed molecular weight by Western-blot: 58 kDa

Preparation

This antibody was obtained by immunization of DKC1 recombinant protein (Accession Number: NM_001363). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109981(DKC1 Antibody) at dilution of 1:600