

# DKC1 antibody



Catalog Number: 109981

## Product name

DKC1 antibody

## Specificity

Human; other species not tested.

## Antibody description

DKC1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, COLO 320 cells, HeLa cells, Jurkat cells, MCF-7 cells. Observed molecular weight by Western-blot: 58 kDa

## Preparation

This antibody was obtained by immunization of DKC1 recombinant protein (Accession Number: NM\_001363). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109981(DKC1 Antibody) at dilution of 1:600