

DLC1 antibody

Catalog Number: 109986

Product name

DLC1 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Immunogen

[Human DLC1 Recombinant protein \(His tag\)](#)

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DLC1 Rabbit Polyclonal antibody. Positive IP detected in mouse lung tissue. Positive IHC detected in human spleen tissue, human liver tissue. Positive IF detected in HepG2 cells.

Preparation

This antibody was obtained by immunization of DLC1 recombinant protein (Accession Number: XM_005273374). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC, IF, IP


Dilutions

Recommended Dilution:


Validations



Immunohistochemical of paraffin-embedded human spleen using Catalog No:109986(DLC1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human spleen using Catalog No:109986(DLC1 antibody) at dilution of 1:100 (under 40x lens)



DLC1 antibody

Catalog Number: 109986

Immunofluorescent analysis of HepG2 cells, using DLC1 antibody Catalog No:109986 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).

IP Result of anti-DLC1 (IP:Catalog No:109986, 4ug; Detection:Catalog No:109986 1:300) with mouse lung tissue lysate 2400ug.

