

DLL4 antibody

Catalog Number: 109997

Product name

DLL4 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

DLL4 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, A431 cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 75-80 kDa

Preparation

This antibody was obtained by immunization of DLL4 recombinant protein (Accession Number: NM_019074). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations




mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109997(DLL4 antibody) at dilution of 1:300



-	+	-	+	Overexpression of DLL4-myc-his
+	+	-	-	Rabbit Anti human DLL4 polyclonal antibody
-	-	+	+	Mouse Anti his-tag monoclonal antibody

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109997(DLL4 Antibody) at dilution of 1:1000



IP Result of anti-DLL4 (IP:Catalog No:109997, 4ug; Detection:Catalog No:109997 1:300) with mouse brain tissue lysate 2800ug.