

DLX1 antibody



Catalog Number: 109998

Product name

Polyclonal

DLX1 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

DLX1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, A375 cells, HeLa cells. Observed molecular weight by Western-blot: 27kd

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of DLX1 recombinant protein (Accession Number: NM_178120). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109998(DLX1 antibody) at dilution of 1:300

