

DMRTC1B antibody



Catalog Number: 110007

Product name

DMRTC1B antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

DMRTC1B Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, mouse spleen tissue. Observed molecular weight by Western-blot: 13 kDa

Preparation

This antibody was obtained by immunization of DMRTC1B recombinant protein (Accession Number: BC029799). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110007(DMRTC1B antibody) at dilution of 1:500