

# Drebrin antibody



Catalog Number: 110019

## Product name

Drebrin antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Drebrin Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, K-562 cells. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 100-130 kDa

## Preparation

This antibody was obtained by immunization of Drebrin recombinant protein (Accession Number: NM\_004395). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:1000-1:10000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110019(DBN1 antibody) at dilution of 1:1500



IP Result of anti-Drebrin (IP:Catalog No:110019, 4ug; Detection:Catalog No:110019 1:2000) with K-562 cells lysate 1200ug.