

Drosha antibody



Catalog Number: 110023

Product name

Drosha antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Drosha Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 90 kDa, 151 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_013235). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110023(Drosha antibody) at dilution of 1:800