

# DSG4 antibody



Catalog Number: 110028

## Product name

DSG4 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

DSG4 Rabbit Polyclonal antibody. Positive WB detected in mouse skin tissue. Observed molecular weight by Western-blot: 160 kDa

## Preparation

This antibody was obtained by immunization of DSG4 recombinant protein (Accession Number: NM\_177986). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:110028(DSG4 Antibody) at dilution of 1:600