

DSG4 antibody



Catalog Number: 110028

Product name

Polyclonal

DSG4 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

DSG4 Rabbit Polyclonal antibody. Positive WB detected in mouse skin tissue. Observed molecular weight by Western-blot: 160 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of DSG4 recombinant protein (Accession Number: NM_177986). Purification method: Antigen Affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:110028(DSG4 Antibody) at dilution of 1:600