

DTD1 antibody

Catalog Number: 110031

Product name

DTD1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

DTD1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse spleen tissue, mouse testis tissue. Observed molecular weight by Western-blot: 23-28kd

Preparation

This antibody was obtained by immunization of DTD1 recombinant protein (Accession Number: NM_080820). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110031(DTD1 antibody) at dilution of 1:500