

# DTNB antibody

Catalog Number: 110034

---

## Product name

DTNB antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

DTNB Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 64kd

## Preparation

This antibody was obtained by immunization of DTNB recombinant protein (Accession Number: NM\_183361). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:110034(DTNB antibody) at dilution of 1:300