

DTX3 antibody



Catalog Number: 110038

Product name

DTX3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

DTX3 Rabbit Polyclonal antibody. Positive WB detected in mouse ovary tissue, mouse brain tissue, mouse kidney tissue, mouse testis tissue. Positive IHC detected in human kidney tissue, human colon tissue. Observed molecular weight by Western-blot: 38 kDa

Preparation

This antibody was obtained by immunization of DTX3 recombinant protein (Accession Number: XM_047428484). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:110038(DTX3 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:110038(DTX3 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog

DTX3 antibody



Catalog Number: 110038

No:110038(DTX3 Antibody) at dilution of 1:50 (under 40x lens)