

# DTX3L antibody



Catalog Number: 110039

## Product name

DTX3L antibody

## Specificity

human; other species not tested.

## Antibody description

DTX3L Rabbit Polyclonal antibody. Positive WB detected in human lung tissue. Observed molecular weight by Western-blot: 84kd

## Preparation

This antibody was obtained by immunization of DTX3L recombinant protein (Accession Number: NM\_138287). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110039(DTX3L antibody) at dilution of 1:500