

DPEP1 antibody

Catalog Number: 110062

Product name

DPEP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DPEP1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue, human kidney tissue. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of DPEP1 recombinant protein (Accession Number: NM_001389466). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:110062(DPEP1 antibody) at dilution of 1:50



(under 10x lens)

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110062(DPEP1 Antibody) at dilution of 1:300