

# DPEP2 antibody

Catalog Number: 110063

## Product name

DPEP2 antibody

## Immunogen

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DPEP2 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 45 kDa, 50 kDa

## Preparation

This antibody was obtained by immunization of DPEP2 recombinant protein (Accession Number: BC024021). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110063(DPEP2 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:110063(DPEP2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:110063(DPEP2 Antibody) at dilution of 1:50

# DPEP2 antibody



Catalog Number: 110063

---

(under 40x lens)