

DPPA2 antibody



Catalog Number: 110073

Product name

DPPA2 antibody

Specificity

Human; other species not tested.

Antibody description

DPPA2 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 34kd

Preparation

This antibody was obtained by immunization of DPPA2 recombinant protein (Accession Number: XM_011512445). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110073(ECSA antibody) at dilution of 1:500