

# DNAJA2 antibody



Catalog Number: 110086

## Product name

Polyclonal

DNAJA2 antibody

## Ig Type

## Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

DNAJA2 Rabbit Polyclonal antibody. Positive WB detected in A375 cells, human brain tissue, PC-3 cells. Observed molecular weight by Western-blot: 46kd

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Preparation

This antibody was obtained by immunization of DNAJA2 recombinant protein (Accession Number: NM\_005880). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:110086(DNAJA2 antibody) at dilution of 1:500