

DULLARD antibody



Catalog Number: 110110

Product name

DULLARD antibody

Specificity

Human; other species not tested.

Antibody description

DULLARD Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of DULLARD recombinant protein (Accession Number: NM_015343). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110110(DULLARD Antibody) at dilution of 1:300