

DUSP22 antibody

Catalog Number: 110119

Product name

DUSP22 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

DUSP22 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells, mouse kidney tissue. Positive IHC detected in human lung cancer tissue. Observed molecular weight by Western-blot: 21 kDa

Preparation

This antibody was obtained by immunization of DUSP22 recombinant protein (Accession Number: NM_020185). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:110119(DUSP22 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:110119(DUSP22 antibody) at dilution of 1:50 (under 40x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110119(DUSP22 Antibody) at dilution of

DUSP22 antibody



Catalog Number: 110119

1:1000