

Dymeclin antibody



Catalog Number: 110128

Product name

Dymeclin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Dymeclin Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 76-80kd

Preparation

This antibody was obtained by immunization of Dymeclin recombinant protein (Accession Number: BC001252). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110128(DYM antibody) at dilution of 1:600