

E2A antibody

Catalog Number: 110152

Product name

E2A antibody

Immunogen

[Human TCF3 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

E2A Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, K-562 cells, PC-3 cells. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 68 kDa

Preparation

This antibody was obtained by immunization of E2A recombinant protein (Accession Number: XM_017027180). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110152(TCF3 antibody) at dilution of 1:500



IP Result of anti-TCF3 (IP:Catalog No:110152, 4ug; Detection:Catalog No:110152 1:500) with Jurkat cells lysate 2800ug.