

E2F7 antibody



Catalog Number: 110156

Product name

E2F7 antibody

Specificity

Human; other species not tested.

Antibody description

E2F7 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IHC detected in human cervical cancer tissue. Observed molecular weight by Western-blot: 95 kDa

Preparation

This antibody was obtained by immunization of E2F7 recombinant protein (Accession Number: NM_203394). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110156(E2F7 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using Catalog No:110156(E2F7 Antibody) at dilution of 1:50 (under 40x lens)