

# E2F8 antibody

Catalog Number: 110157

## Product name

E2F8 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

E2F8 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 82 kDa, 105 kDa

## Preparation

This antibody was obtained by immunization of E2F8 recombinant protein (Accession Number: BC028244). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110157(E2F8 antibody) at dilution of 1:1000



IP Result of anti-E2F8 (IP:Catalog No:110157, 4ug; Detection:Catalog No:110157 1:500) with HEK-293 cells lysate 2800ug.