

EAF1 antibody

Catalog Number: 110161

Product name

EAF1 antibody

IP: 1:200-1:2000

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EAF1 Rabbit Polyclonal antibody. Positive IF detected in SH-SY5Y cells. Positive IP detected in HeLa cells. Positive WB detected in K-562 cells, HeLa cells, mouse brain tissue. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of EAF1 recombinant protein (Accession Number: NM_033083). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IP, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110161(EAF1 antibody) at dilution of 1:500



IP Result of anti-EAF1 (IP:Catalog No:110161, 4ug; Detection:Catalog No:110161 1:500) with HeLa cells lysate 1200ug.



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:110161(EAF1 Antibody) at

EAF1 antibody



Catalog Number: 110161

dilution of 1:50 and Alexa Fluor 488-conjugated

AffiniPure Goat Anti-Rabbit IgG(H+L)