

EAF6 antibody

Catalog Number: 110163

Product name

EAF6 antibody

Specificity

Human; other species not tested.

Antibody description

EAF6 Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in HepG2 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 22-28 kDa

Preparation

This antibody was obtained by immunization of EAF6 recombinant protein (Accession Number: NM_001270875). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110163(C1orf149 Antibody) at dilution of 1:600



IP Result of anti-C1orf149 (IP:Catalog No:110163, 4ug; Detection:Catalog No:110163 1:500) with HepG2 cells lysate 3480ug.