

# EIF1B antibody



Catalog Number: 110177

## Product name

EIF1B antibody

## Specificity

Human; other species not tested.

## Antibody description

EIF1B Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HeLa cells, HepG2 cells. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 14 kDa

## Preparation

This antibody was obtained by immunization of EIF1B recombinant protein (Accession Number: NM\_005875). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:110177(EIF1B antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:110177(EIF1B Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog

# EIF1B antibody



Catalog Number: 110177

---

No:110177(EIF1B Antibody) at dilution of 1:50 (under 40x lens)