

EIF2B1 antibody

Catalog Number: 110181

Product name

EIF2B1 antibody

IP: 1:200-1:2000

IF: 1:10-1:100

Specificity

Human; other species not tested.

Antibody description

EIF2B1 Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in MCF7 cells, HeLa cells, K-562 cells. Positive IF detected in Hela cells. Observed molecular weight by Western-blot: 34 kDa

Preparation

This antibody was obtained by immunization of EIF2B1 recombinant protein (Accession Number: NM_001414). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110181(EIF2B1 antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells, using EIF2B1 antibody Catalog No:110181 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent



EIF2B1 antibody



Catalog Number: 110181

IP Result of anti-EIF2B1 (IP:Catalog No:110181,

4ug; Detection:Catalog No:110181 1:700) with
K-562 cells lysate 3200ug.