

# EIF2C3-Specific antibody



Catalog Number: 110186

## Product name

Polyclonal

EIF2C3-Specific antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

EIF2C3-Specific Rabbit Polyclonal antibody.  
Positive WB detected in SKOV-3 cells, HepG2 cells. Observed molecular weight by Western-blot: 90-100 kDa

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_024852).  
Purification method: Antigen affinity purified.

## Formulation


PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

## Validations



SKOV-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110186(EIF2C3-Specific Antibody) at dilution of 1:300