

EIF3A antibody



Catalog Number: 110190

Product name

EIF3A antibody

Specificity

Human; other species not tested.

Antibody description

EIF3A Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 166 kDa

Preparation

This antibody was obtained by immunization of EIF3A recombinant protein (Accession Number: NM_003750). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:50-1:500

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110190(EIF3A Antibody) at dilution of 1:600



Immunofluorescent analysis of (10% Formaldehyde) fixed HEK-293 cells using Catalog No:110190(EIF3A Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)