

EIF3J antibody

Catalog Number: 110198

Product name

EIF3J antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

IF: 1:20-1:200

Antibody description

EIF3J Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A549 cells, HEK-293 cells, Jurkat cells, K-562 cells, MCF7 cells, SKOV-3 cells. Positive IP detected in HeLa cells. Positive IF detected in HeLa cells. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 35 kDa, 60-66 kDa

Validations



Preparation

This antibody was obtained by immunization of EIF3J recombinant protein (Accession Number: NM_003758). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110198(EIF3J antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IF, IP

Dilutions


Recommended Dilution:



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:110198(EIF3J antibody) at dilution of 1:100 (under 10x lens)


EIF3J antibody

Catalog Number: 110198



Immunofluorescent analysis of HeLa cells, using EIF3J antibody Catalog No:110198 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).

Blue pseudocolor = DAPI (fluorescent DNA dye).



IP Result of anti-EIF3J (IP:Catalog No:110198, 4ug; Detection:Catalog No:110198 1:500) with HeLa cells lysate 880ug.