

EFCAB1 antibody

Catalog Number: 110203

Product name

EFCAB1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EFCAB1 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of EFCAB1 recombinant protein (Accession Number: NM_024593). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110203(EFCAB1 Antibody) at dilution of 1:1000