

EFCAB1 antibody



Catalog Number: 110203

Product name

Polyclonal

EFCAB1 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

EFCAB1 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 30 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of EFCAB1 recombinant protein (Accession Number: NM_024593). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110203(EFCAB1 Antibody) at dilution of 1:1000