

EFCAB4B antibody



Catalog Number: 110204

Product name

EFCAB4B antibody

Specificity

Human; other species not tested.

Antibody description

EFCAB4B Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 46kd

Preparation

This antibody was obtained by immunization of EFCAB4B recombinant protein (Accession Number: NM_032680). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110204(EFCAB4B Antibody) at dilution of 1:600