

EFCAB5 antibody

Catalog Number: 110205

Product name

EFCAB5 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human; other species not tested.

Validations

Antibody description

EFCAB5 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, L02 cells. Positive IHC detected in human liver cancer tissue, human colon tissue, human testis tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 110 kDa




Preparation

This antibody was obtained by immunization of EFCAB5 recombinant protein (Accession Number: NM_198529). Purification method: Antigen affinity purified.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110205(EFCAB5 Antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human liver cancer slide using Catalog No:110205(EFCAB5 Antibody) at dilution of 1:50

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

Immunohistochemistry of paraffin-embedded human liver cancer slide using Catalog

EFCAB5 antibody

Catalog Number: 110205

No:110205(EFCAB5 Antibody) at dilution of 1:50



Immunofluorescent analysis of HepG2 cells using Catalog No:110205(EFCAB5 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG