

EFCAB5 antibody



Catalog Number: 110205

Product name

EFCAB5 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human; other species not tested.

Antibody description

EFCAB5 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, L02 cells. Positive IHC detected in human liver cancer tissue, human colon tissue, human testis tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 110 kDa

Preparation

This antibody was obtained by immunization of EFCAB5 recombinant protein (Accession Number: NM_198529). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110205(EFCAB5 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human liver cancer slide using Catalog No:110205(EFCAB5 Antibody) at dilution of 1:50




Immunohistochemistry of paraffin-embedded human liver cancer slide using Catalog

EFCAB5 antibody



Catalog Number: 110205

No:110205(EFCAB5 Antibody) at dilution of 1:50



Immunofluorescent analysis of HepG2 cells using Catalog No:110205(EFCAB5 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG