

# EMG1 antibody



Catalog Number: 110238

## Product name

EMG1 antibody

IP: 1:200-1:1000

IF: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EMG1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, HEK-293 cells, PC-3 cells. Observed molecular weight by Western-blot: 27 kDa

## Preparation

This antibody was obtained by immunization of EMG1 recombinant protein (Accession Number: NM\_006331). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IF, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110238(EMG1 antibody) at dilution of 1:400



IP Result of anti-EMG1 (IP:Catalog No:110238, 4ug; Detection:Catalog No:110238 1:400) with HeLa cells lysate 1080ug.



Immunofluorescent analysis of HeLa cells using Catalog No:110238(EMG1 Antibody) at dilution of

# EMG1 antibody



Catalog Number: 110238

---

1:50 and Alexa Fluor 488-conjugated AffiniPure

Goat Anti-Rabbit IgG(H+L)