

EML1 antibody

Catalog Number: 110241

Product name

EML1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

EML1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse brain tissue. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 90-92kd

Preparation

This antibody was obtained by immunization of EML1 recombinant protein (Accession Number: NM_004434). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110241(EML1 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using EML1 antibody Catalog No:110241 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).