

ENC1 antibody

Catalog Number: 110246

Product name

ENC1 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human,Mouse,Rat; other species not tested.

Antibody description

ENC1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas tissue. Positive WB detected in mouse brain tissue, human lung tissue, MCF7 cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 67 kDa,57 kDa

Preparation

This antibody was obtained by immunization of ENC1 recombinant protein (Accession Number: NM_001256575). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110246(ENC1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:110246(ENC1 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human pancreas using Catalog No:110246(ENC1

ENC1 antibody

Catalog Number: 110246

antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-ENC1 (IP:Catalog No:110246, 4ug; Detection:Catalog No:110246 1:1000) with mouse brain tissue lysate 2640ug.