

ENO3 antibody

Catalog Number: 110255

Product name

ENO3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ENO3 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, human brain tissue. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of ENO3 recombinant protein (Accession Number: NM_053013). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:110255(ENO3 antibody) at dilution of 1:500